

Distal Femur Outcomes in High Risk Patients*

Stainless Steel Plate

vs.

Carbon Fiber PEEK Plate




PATIENT OUTCOME
64% BONE UNION

PATIENT OUTCOME
91% BONE UNION




COMPLICATIONS
18% HARDWARE FAILURES

COMPLICATIONS
0% HARDWARE FAILURES



PATIENT RISK FACTORS
9% DIABETES
27% SMOKER
0% PVD+
57 AVG. AGE


PATIENT RISK FACTORS
55% DIABETES
18% SMOKER
18% PVD+
71 AVG. AGE



CLINICAL EXPERIENCE

40 Publications covering 731 Cases*

ANATOMIC LOCATIONS AND NUMBER OF CASES



*As of January 2020.

REFERENCES

- Hak DJ, et al. Paper 482 presented at: EFORT 2014 Annual Meeting; 2014 Jun 4-6; London, UK
- Hak DJ, et al. *Injury*. 2018 Jun;49 Suppl 1:S8-S11.
- Krantzow M. *JOT Case Report* 2019 <https://journals.lww.com/jorthotrauma/Documents/JOT-2019-324.pdf>
- Gurich RW, et al. *Orthopedics Today* 2018 Oct. <https://www.healio.com/orthopedics/oncology/news/print/orthopedics-today/%7Be20443bf-31c1-420d-87cb-ceb91dc92777%7D/carbon-fiber-implants-aid-fixation-of-extremity-fractures-in-oncology-cases>
- Hillock R, et al. *JISRF Reconstructive Review*. 2014 Mar;4(1):23-32.
- Rosenthal HG. Paper 35 presented at: MSTs 2018 Annual Meeting; 2018 Oct 10-12; New York, NY.
- Hak DJ, et al. *Orthopedics*. 2014 Dec;37(12):825-30.
- Sacchetti F, et al. *Eur J Orthop Surg Traumatol*. 2019 Dec 21. <https://link.springer.com/article/10.1007%2Fs00590-019-02602-4>
- Piccoli A, et al. *Injury*. 2017 Oct;48 Suppl3:S55-S59.
- Vles GF, et al. *Case Rep Orthop*. 2019 Jul 30; 2019:1538158.
- Laux CJ, et al. *World J Surg Oncol*. 2018 Dec 28;16(1):241.
- Collis PN, et al. *Injury*. 2011 Apr;42(4):424-6.
- Caforio M, et al. *J J Ortho Rheum*. 2015, 1(2): 008.
- Kojic N, et al. *Injury*. 2017 Nov;48 Suppl 5:S8-S11.
- Rives AW. Case Report. <https://www.carbo-fix.com/publications-czs>
- Merolli A, et al. *J Mater Sci Mater Med*. 2016 Mar;27(3):50.
- Di Maggio B, et al. *Injury*. 2017 Oct;48 Suppl 3:S34-S38.
- Perugia D, et al. *Injury*. 2017 Oct;48 Suppl 3:S24-S29.
- Guzzini M, et al. *Acta Biomed*. 2018 Oct 23;90(1-5):61-66.
- Zoccali C, et al. *J Orthop*. 2016 Jul 6;13(4):343-6.
- Fischler T. *JOT Case Report* 2019. <https://journals.lww.com/jorthotrauma/Documents/JOT10246.pdf>
- Prince DE, et al. Poster presented at: MSTs 2013 Annual Meeting; 2013 Oct 3-5; San Francisco, CA. <https://www.carbo-fix.com/publications-czs>
- Ziran BH, et al. Poster 68 presented at: OTA 2017 Annual Meeting; 2017 Oct 11-14; Vancouver, Canada.
- Ziran BH, et al. *JOT Case Report* 2018. <https://journals.lww.com/jorthotrauma/Documents/JOT9865.pdf>
- Vercio RC, et al. *J Am Acad Orthop Surg*. 2019 Jun 15;27(12):e585-e588.
- Allison DC, et al. Poster presented at: MSTs 2013 Annual Meeting; 2013 Oct 3-5; San Francisco, CA. <https://www.carbo-fix.com/publications-czs>
- Zimel MN, et al. *Skeletal Radiol*. 2015 Sep;44(9):1317-25.
- Fragomen AT, et al. *HSS J*. 2019 Jul;15(2):109-114.
- Lee AK, et al. Paper 155 presented at: OTA 2017 Annual Meeting; 2017 Oct 11-14; Vancouver, Canada.
- Fragomen AT, et al. *HSS J*. 2018 Oct;14(3):258-265.
- Niemann J. Case Study Video. <https://www.youtube.com/watch?v=3kNuCax3ytM>
- Ortega GR et al. *JOT Case Report* 2018. <https://journals.lww.com/jorthotrauma/Documents/JOT9970.pdf>
- Loeb AE, et al. *JBJS Case Connect*. 2018 Oct-Dec;8(4):e83.
- Ziran BH, et al. Carbon Fiber Reinforced PEEK Versus Titanium Tibial Intramedullary Nailing: A Preliminary Analysis and Results. In submission. *JOT*.
- Ortega GR. Digital Case Presentation 15 presented at: OTA 2019 Annual Meeting; 2019 Sep 25-28. Denver, CO.
- Caforio M, et al. *Injury*. 2014 Dec;45 Suppl 6:S36-8.
- Guzzini M, et al. *Injury*. 2017 Jun;48(6): 1249-1252.
- Pinter ZW, et al. *Foot Ankle Spec*. 2018 Jun;11(3):223-229.
- Canales MB, et al. *Podiatry Today* 2017 June;30(6). <https://www.podiatrytoday.com/exploring-carbon-fiber-fixation-lower-extremity>
- Araoye IB, et al. *Foot Ankle Surg*. 2019 Feb;25(1):2-7.

For further information call us toll free at 866-INVIBIO or +44 (0)1253 898000 or please visit our website at:

▶ **Invibio.com**

Victrex plc and/or its group companies (Victrex plc) believes that the information contained in this document is an accurate description of the typical characteristics and/or uses of the product(s) and is based on information that we believe is reliable. However, it is provided for information only. It is not intended to amount to advice on which you should rely and should not be construed as, or used as a substitute for, professional medical advice or other professional or specialist advice. In particular, it is the customer's responsibility to thoroughly test the product in each specific application to determine its performance, efficacy, and safety for each end-use product, device or other application. Suggestions of product uses should not be taken as inducements to infringe any particular patent. Mention of a product in this document is not a guarantee of its availability.

Victrex plc reserves the right to modify products, specifications and/or packaging as part of a continuous program of product development. Victrex plc makes no warranties, express or implied, including, without limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, in fact or by law.

Further, Victrex plc makes no warranty to your customers or agents and has not authorized anyone to make any representation or warranty other than as provided above. Victrex plc shall in no event be liable for any general, indirect, special, consequential, punitive, incidental or similar damages, or any damages for harm to business, lost profits or lost savings, even if Victrex has been advised of the possibility of such damages regardless of the form of action. The foregoing does not seek to affect any liability (including to individual consumers) which cannot be excluded or limited under any applicable law.

Supporting information is available on request for claims referenced in this document.

Victrex plc is the owner or the licensee of all intellectual property rights in this document. All rights are protected by copyright and other intellectual property laws and treaties around the world. All rights reserved.

VICTREX™ and SHAPING FUTURE PERFORMANCE™ are trademarks of Victrex plc or its group companies.

INVIBIO™, JUVORA™, PEEK-OPTIMA™, BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies.